In The Name Of God

Shaheed Beheshti University of
Medical sciences and Health Services

Collaborating Center for Tuberculosis
East Mediterranean Region

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National Research Institute
Tuberculosis and Lung Disease

Collaborating Center for Tuberculosis
Laboratory IUATLD

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Introduction

The Iranian Reference TB Laboratory receives between 80 to 120 pulmonary and Extra-Pulmonary specimens per day. The specimens were digested and decontaminated by petroff's (using 4% NAOH) method.

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Specimen's collection

Specimens	Collection instruction (time & Preservatives)	Volume/ml	Specimens storage in 4^0C	Unacceptable Specimens
Abscess contents, aspirate fluid	Clean the skin with alcohol before aspirating .Swab can placed in small amount of 7H9 medium for transport	-	8-12 hours	Dry swab
Blood	SPS (sodium polyanethol sulfonate) or Heparin as anticoagulant	5-10 ml	8-12 hours	-
Body fluids (pleural, pericardial ,peritoneal, etc)	Bloody specimens collected in 1% Heparin	10-15 ml	Maximum 24 hours	EDTA inhibits mycobacterial growth even in trace amounts
Bone marrow	Bloody specimens collected in 1% Heparin	1-2 ml	Maximum 24 hours	EDTA inhibits mycobacterial growth even in trace amounts
Bronchoalveolar lavage or bronchial washings	Collect without preservative	5-10 m	8-12 hours	EDTA inhibits mycobacterial growth even in trace amounts
CSF	Collect without preservative	2-5 ml	8-12 hours	
Gastric lavage fluid	Collect fasting early-morning specimens on 3 consecutive days. Adjust the PH with (100 mg/10 ml specimens) of sodium carbonate immediately after collecting.	5-10 ml	2-4 hours	
Lymph node	Collect without preservative.	-	Maximum 24 hours	don't immerse in saline or other fluid
Sputum (expectorated)	Rinse the mount and collect 3 consecutive morning sputum	5-10 ml	3-5 days	
Sputum (induced)	Use sterile hypertonic saline	-	8-12 hours	
Stool	Collect directly into container	-	Maximum 24 hours	
Tissue biopsy sample/wound material	Collect without preservative	1g	Maximum 24 hours	don't immerse in saline or other fluid
Urine	Collect first morning specimen on 3 consecutive days	minimum 40 ml		24-h pooled specimens; urine from catheter bag

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Specimens should be transported to the laboratory to avoid overgrowth by other microorganisms. This is particularly true for specimens from non-sterile-sites such as sputum. When the transport or the processing is delayed, specimens should be stored not more than 5 days at 4^0C. The cetylpyridinium chloride (CPC) method is widely used for the transport of sputum specimens .CPC eliminates the associated flora in sputum specimens and treated specimens should not be submitted to further decontamination prior to cultivation. Specimens treated with CPC preferentially inoculated in egg-based culture media.

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Cetylpyridinium Chloride (CPC) a quaternary ammonium compound is used to decontaminate the specimens while sodium chloride effects liquefaction.

Materials

1% cetylpyridinium chloride, 2% sodium chloride, D.W or sterile saline.

Preparation

Dissolve 10 g of cetylpyridinium chloride and 20 g of Nacl in 1000 ml D.W .Solution should be protected from light, heat and evaporation.

Procedures:

To 10 ml-of sputum specimens, add equal volume of CPC-Nacl solution .Shake by hand. Package specimens in mailing container and send to the laboratory for processing.

Once received in the laboratory, dilute the specimens with D.W. to 50ml- mark. Centrifuges at 3000 x g for 15 minutes Decant supernatant fluid and suspend sediment in 1 to 2 ml of D.W, saline or 0.2% bovine albumin fraction. The re suspended medium inoculated onto egg –medium and incubated at 35^0 to 37^0C.

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The rapidly growing normal flora contaminate majority of clinical specimens that submitted to the tuberculosis culture laboratory. They would rapidly overgrow the entire surface of the medium and digest it before the tubercle bacilli start to grow. Therefore most specimens must be subjected to a harsh digestion and decontamination procedures that liquefies the organic debris and eliminates the unwanted flora.

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When culturing tubercle bacilli, two important aspects should be borne in mind.
Specimens homogenized to free the bacilli from mucus, cells or tissue in which they embedded.
Neither homogenization nor decontamination should unnecessarily diminish the viability of tubercle bacilli.

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Materials for sodium hydroxide (Modified Petroff) method

4% NaOH solution, 2 NHCl solutions –phenol red indicators

Preparation

NaOH: Dissolve 4 g NaOH in 100ml distilled water and sterilize by autoclaving at 121^0C for 15 minutes (4% NaOH)

HCL: Prepare the 2.0 N HCL by carefully adding 33ml concentrated HCL to 200 ml distilled water and sterilize by autoclaving


Phenol-red indicator: 80 mg phenol red powder in 20 ml 4% NaOH and add sufficient distilled water to make 1000 ml

HCL- Phenol red indicator: combine 20 ml of phenol red solution (80 mg phenol red powder in 20 ml 4% NaOH) and 85 ml of concentrated HCL, with distilled water to make 1000ml

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To equal volume of sputum specimen add equal amount of 4% NaOH.

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Tighten the container cap and shake to digest.

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Let stand for 15 minutes at room temperature.

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
Neutralize sediment* with HCl until indicator changes from red to yellow. Re-suspend the sediment with 1 to 2 ml of sterile water.

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Mix sediment and inoculate it to L.J culture medium, (For extra-pulmonary specimen three L.J media and for sputum two L.J media are used)

*. Adjusting the PH of the sediment is not necessary because of the well –buffered media. But do not over-inoculate because this interferes with the buffer action and results in poor growth.

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Materials


5% Oxalic –acid
4% Sodium hydroxide
Physiological saline
Indicator solution, phenol red

Preparation

1) Oxalic acid ($C_2H_2O_4.2H_2O$): mix 5 g per 100 ml D.W and autoclave.
2) Sodium hydroxide: Mix 4 g per 100 ml distilled water and autoclave.
3) Physiological saline: Mix 0.85 g per 100 ml distilled water and autoclave.
4) Phenol red indicator: Dissolve 80 mg phenol red powder in 20 ml 4% NaOH, bring volume up to 1,000 ml with D.W.

Procedure

1) Add an equal volume of 5% oxalic acid to each specimen.
2) Mix on test –tube mixer. Let stand for 30 minutes with occasional shaking.
3) Add sterile saline to the 50-ml mark on the centrifuge tube.
4) Centrifuge the specimen at 2,000 to 3,000 X g for 15 minutes.
5) Decant the supernatant fluid into splash-proof container.
6) Add a few drops of the phenol-red indicator to the sediment and neutralize with the 4% NaOH until the first persistent pale pink color forms.
7) Inoculate media and incubate at 35^0C to 37^0C.

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
LÖWENSTEIN-JENSEN MEDIUM

1) Ingredients: Mineral salt solution, Malachite green solution 2%, homogenised whole eggs.

2) Preparation of complete medium

Coagulation of medium, Sterility check, Storage

Mineral salt solution

Potassium dihydrogen phosphate anhydrous	2.4g
Magnesium sulphate	0.24g
Magnesium citrate	0.6g
Asparagine	3.6g
Glycerol	12 ml
Distilled water	600 ml

Dissolve the ingredients in order in the distilled water by heating. Autoclave at 121 ºC for 30 minutes to sterilize.

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Malachite green solution 2%:

Malachite green dye	2 g
Sterile distilled water	100 ml

Dissolve the dye in sterile distilled water by placing the solution on stirrer. Distribute the solutions in special tubes, then autoclave it. Always keep the malachite solution in dark.

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Homogenised whole eggs

Clean fresh eggs by scrubbing thoroughly with a hand brush in warm and a plain alkaline soap. Rinse eggs in running water and soak them in 70% ethanol for 15 minutes. Before handling the clean dry eggs scrub the hands and wash them. Crack the eggs with a sterile knife into a sterile flask and beat them with a sterile egg whisk or in a sterile blender.

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Preparation of complete medium

The following ingredients are aseptically pooled in a large, sterile flask and mixed well:

Mineral salt solution:	600 ml
Malachite green solution:	20 ml
Homogenised eggs:	1000 ml

Complete egg medium is distributed in 6-8 volumes in sterile McCartney bottles and the tops are tightened.

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Coagulation of medium

Before loading, heat the inspissator to 80 ºC to quicken the build-up of the temperature. Place the bottles in a slanted position in the inspissator and coagulate the medium for 45 minutes at 80º-85º C (Since the medium has been prepared with sterile precautions this heating is to solidify the medium, not to sterilise it). Heating for a second or third time has a detrimental effect on the quality of the medium.

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Middle brook 7 H-10 Agar Medium

Middle brook 7H-10 agar medium may be prepared by using commercially available 7H-10 agar – powdered base and Middle brook OADC enrichment.

To prepare 200 ml of complete medium

Suspend 3.6 g of Middlebrook 7H-10 basal medium in 180 ml of sD.W snd add 1.0 ml of glycerol (reagent grade).

Autoclave for 10 minutes at 121^0C. As soon as cooled to 50 to 56^oC, add 20 ml of OADC. Dispense within 1 hour after autoclaving. Allow to solidify at room temperature without exposure to daylight.

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Middlebrook 7 H-9 Liquid Medium

Middle brook 7H-10 agar medium may be prepared by using commercially available 7H-10 agar – powdered base and Middle brook OADC enrichment.

To prepare 1000 ml of complete medium
Suspend 4.7 g of dehydrated medium in 900 ml of sD.W containing 2 ml of glycerol.
Autoclave for 15 minutes at 121^0 C. As soon as cooled to 45^0 C, add 100 ml of OADC. Dispense 5-ml amounts into 20X 150 mm screw cap tubes.

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Inoculation

Glass or plastic Pasture pipettes used. Each slope inoculated with 0.2-0.4 ml (2-4 drops) of the centrifuged sediment.

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Incubation

All cultures incubated at 35º-37º C until growth observed or discarded as negative after eight weeks. Inoculated media preferably incubated in a slanted position for at least 24 hours to ensure the complete distribution of inoculum.

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Examination schedule

Examine all cultures 72 hours after inoculation to check that liquid has completely evaporated.Thereafter, the caps


Tighten in order to prevent drying out of media and to detect contaminants.

- After one week to detect rapidly growing mycobacteria
- Weekly through eight weeks to detect positive cultures of M.tuberculosis

Reading of cultures

Typical colonies of *M.tuberculosis* are rough, crumbly, waxy, non-pigmented (Cream colored) and slow-growers, i.e. only appearing three weeks after inoculation. With doubtful cases, the acid-fastness should be confirmed by Ziehl-Neelsen staining.

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Reading Reporting

No growth	Negative
<50 colonies	Number of colonies
50-100 colonies	Positive (+1)
100-200 colonies	Positive (+2)
200-500 colonies	Positive (+3)
> 500 colonies	Positive (+4)
Contamination	Contamination

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The definitive diagnosis of mycobacterial disease demands that the causative agent recovered on culture medium.
Culture techniques have been estimated to detect as many as 10 to 1000 viable mycobacteria per ml. of specimens. A positive culture for *M.tuberculosis* confirms the diagnosis of active disease.

Address: Mycobacteriology Research Centre, Iranian National Reference TB Laboratory, National Research Institute of Tuberculosis & Lung Disease (NRITLD), Shaheed Beheshti University of Medical Sciences, Shaheed Bahonar Ave, Darabad, 19556, P.O 19575/154, Tehran, Iran.
Tel: (0098-21) 20105555 or 20109505
Fax: (0098-21) 20109484 or 20109505
E-mail: mrc @ nritld.ac.ir, pfarnia @ hotmail.com, pfarnia @ nritld.ac.ir, Website: www.mrc.ac.ir

This document was created with Win2PDF available at http://www.daneprairie.com.
The unregistered version of Win2PDF is for evaluation or non-commercial use only.