In The Name Of God

Shaheed Beheshti University of
Medical sciences and Health Services

Collaborating Center for Tuberculosis
East Mediterranean Region

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National Research Institute
Tuberculosis and Lung Disease

Collaborating Center for Tuberculosis
Laboratory IUATLD

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Drug Susceptibility Test against Mycobacterium tuberculosis

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Introduction

The Drug Susceptibility Test was performed by L.J proportional method on positive cultures.

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The proportional method enables precise quantization of the proportion of mutants resistant to a given drug. Several dilutions of inoculums are planted in replicate onto both control and drug containing media. The proportion of bacilli resistant to given drug is then determined by comparing these numbers and expressing the resistant portion as a percentage of the total population tested. In every case the inoculums size must be known so the actual proportion of resistant mutants may be calculated.

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
Ptoency

The true potency of the drug is the number of micrograms of active drug per milligram total weight of the product. Not all antimicrobial drugs have been isolated in pure form, and a portion of their weight may be due to impurities or to the sulfate or other radical component of the molecule.
To calculate the weight of drug necessary to prepare 10 ml of a 10,000 µg/ml solution, use this formula:

$$\frac{\textbf{10000}}{\textbf{Potency}} \times \textbf{10} = \textbf{milligrams to weight}$$

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Preparation of Isonizid solution


* Dissolve 50 mg of isonizid in 25 ml of distilled water.
 (INH concentration is 2000 mcg /ml = stock solution)
* Add 1 ml of stock solution to 9 ml of distilled water.
 (INH concentration is 200 mcg /ml=solution A)
* Add 1 ml of solution A to 9 ml of distilled water
 (INH concentration is 20 mcg/ml =solution B)
Sterilize solution B by using 0.22 µl Millipore - filter.
* Add 1 ml of filtrated solution B, to 99 ml of L.J culture Media.
 (INH concentration is 0.2 mcg/ml = working solution)

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Preparation of Streptomycin

Preparation of Streptomycin Solution (10 mcg/ml)

Potency of streptomycin sulfate is may be different from company to company.Therefore, each time the weight of drug has to be calculated using this formula:

$$\text{Weight (mg)} = \frac{\text{Volume (ml)X Concentration}}{\text{Ptotency (µg/ml)}}$$

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$$128.7\ (\text{mg}) = \frac{10\times 10{,}000}{777}$$

Preparation of streptomycin solution

* Dissolve 128.7 mg of isonizid in 10 ml of distilled water.

(SM concentration is 10000 mcg /ml= stock solution)

* Add 1 ml of stock solution to 9 ml of distilled water (SM)

(SM concentration is 1000 mcg /ml= solution A)

* Sterilize solution A by using 0.22 µl Millipore -filter
* Add 1 ml of solution A to 996 ml of L.J culture media

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
Preparation of ethambutol


Preparation of ethambutol solution (2mcg/ml)

* Dissolve 50 mg of ethambutol in 25 ml of distilled water.
 (ETB concentration is 2000 mcg /ml= stock solution)
* Add 1 ml of stock solution to 9 ml of distilled water.
 (ETB concentration is 200 mcg /ml= solution A)
* Sterilize solution C by using 0.22 µl Millipore –filter.
* Add 1ml of filtrated solution A, to 99 ml of L.J culture Media.
 (ETB concentration is 2 mcg/ml= working solution)

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Preparation of rifampin solution (40mcg/ml)


* Dissolve 40 mg of rifampin in 10 ml of DMSO or 0.1N HCL
 (RF concentration is 4000 mcg /ml= stock solution)
* Sterilize stock-solution by using 0.22 µl Millipore -filter
* Add 1ml of filtrated solution A, to 99 ml of L.J culture Media
 (RF concentration is 40 mcg/ml= working solution)

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Preparation of Pyrazinamide

Culture Medium preparation:

The commercially available 7H10 agar bases was dissolved in deionized water at 9.5 g per 450 ml .Then 3.3 g of monopotassium phosphate (KH2PO4) was added to acidify the medium. In addition 0.5 g of casein hydrolysate and 2.5 ml of glycerol were added. After being autoclave at 121^0C for 10 min, the medium cool down to 54^0C and 50 ml of OADC added into it. The medium split into to sterile flask and in one of the flask 2.5 ml of pyrazinamide stock solution was added into it.

Preparation of Pyrazinamide solution

* Dissolve 100 mg of Pyrazinamide in 10 ml of distilled Water.
 (PZA concentration is 10000 mcg /ml = Stock solution)
* Sterilize stock-solution by using 0.22 µl Millipore -filter
* Add 2.5 ml of filtrated solution, to 247.5 ml of 7H10-agar culture Media
 (PAZ concentration is 100 mcg/ml =Working solution)

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
Preparation of Capreomycin

Preparation of Capreomycin solution (10mcg/ml)


* Dissolve 100 mg of capreomycin in 10 ml of distilled water
 (CAP concentration is 10000 mcg /ml= stock solution)
* Add 1 ml of stock –solution into 9 ml of distilled–water
 (CAP concentration is 1000 mcg/ml= solution A)
* Sterilize stock-solution by using 0.22 µl Millipore -filter
* Add 1ml of filtrated solution, to 99 ml of L.J culture Media
 (CAP concentration is 10 mcg/ml= working solution)

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Preparation of Ciprofloxacin


Preparation of Ciprofloxacin solution (2mcg/ml)

∗ Dissolve 50 mg of ciprofloxacin in 25 ml of distilled water
(CIP concentration is 2000 mcg /ml = stock solution)
∗ Add 1 ml of stock –solution into 9 ml of distilled –water
(CIP concentration is 200 mcg/ml, solution A)
∗ Sterilize stock-solution by using 0.22 µl Millipore -filter
∗ Add 1ml of filtrated, to 99 ml of L.J culture Media
(CIP concentration is 2 mcg/ml = working solution)

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Preparation of cycloserine solution (30mcg/ml)


* Dissolve 30 mg of cycloserine in 10 ml of distilled water with PH=10 (cycloserine concentration is 3000 mcg/ml= cycloserine)
* Sterilize stock-solution by using 0.22 µl Millipore -filter
* Add 1ml of filtrated solution, to 99 ml of L.J culture Media
(cycloserine concentration is 30 mcg/ml= working solution)

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Preparation of Ethionamide

Preparation of Ethionamide Hydrochloride solution (20mcg/ml)

Potency of ethionamide Hydrochloride is may be different from company to company. Therefore, each time weight of drug necessary has to be calculated using this formula

$$\text{Weight (mg)} = \frac{\text{Volume (ml)} \times \text{Concentration}}{\text{Ptotency (µg/ml)}}$$

To make a 10 ml of 10,000mcg /ml stock solution then

$$= \frac{10{,}000 \times 10}{\text{Potency (mg)}}$$

Dissolve 100 mg of in 10 ml of distilled water
(ETH concentration is 10000 mcg /ml= stock solution)


* Add 1 ml of stock –solution into 9 ml of distilled –water
(ETH concentration is 1000 mcg/ml = solution A)
* Sterilize stock-solution by using 0.22 µl Millipore -filter
* Add 2 ml of filtrated solution A, to 98 ml of L.J culture Media
(ETH concentration is 20 mcg/ml = working solution)

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Preparation of Kanamycin

Preparation of Kanamycin Solution (20 mcg/ml)

* Potency of Kanamycin sulfate is may be different from company to company. Therefore, each time weight of drug necessary has to be calculated using this formula

$$\text{Weight(mg)} = \frac{\text{Volume(ml)} \times \text{Concentration}}{\text{Ptotency (µg/ml)}}$$

To make a 10 ml of 10,000mcg /ml stock solution then $\frac{10{,}000 \times 10}{786} = 127.2$

Preparation of kanamycin solution

* Dissolve 127.2 mg of Kanamycin sulfate in 10 ml of distilled water.
(KM concentration is 10000 mcg /ml= stock solution)
* Add 1 ml of stock solution to 9 ml of distilled water
(KM concentration is 1000 mcg /ml= in solution A)
* Sterilize solution A by using 0.22 µl Millipore -filter
* Add 2 ml of filtrated solution A, to 98 ml of L.J culture Media.

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
Preparation of Ofloxacin

Preparation of Ofloxacin solution (2 mcg/ml)

∗ Dissolve 20 mg of ofloxacin in 10 ml of 0.1NNaOH
(ofloxacin concentration is 20000 mcg/ml= stock solution)

∗ Add 1 ml of stock –solution into 9 ml of distilled –water
(ofloxacin concentration is 2000 mcg/ml= solution A)

∗ Sterilize stock-solution by using 0.22 µl Millipore -filter

∗ Add 1ml of filtrated solution A, to 99 ml of L.J culture Media
(ofloxacin concentration is 2 mcg/ml= working solution)

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
Preparation of P-Amino-salicylic cid

Preparation of P-Amino-salicylic cid Solution (5mcg/ml)

*Potency of P-Amino-salicylic cid is may be different from company to company. Therefore, each time weight of drug necessary has to be calculated. Using this formula

$$\text{Weight (mg)} = \frac{\text{Volume (ml)} \times \text{Concentration}}{\text{Ptotency (µg/ml)}}$$

To make a 10 ml of 10,000mcg /ml stock solution then

$$\frac{10{,}000 \times 10}{\text{(Ptotency) } 877.2}$$

Preparation of kanamycin solution

* Dissolve 114 mg of P-Amino-salicylic acid in 10 ml of distilled water.
 (PAS concentration is 10000 mcg /ml, stock solution)
* Add 1 ml of stock solution to 9 ml of distilled water
 (PAS concentration is 1000 mcg /ml, solution A)
* Sterilize solution A by using 0.22 µl Millipore -filter
* Add 0.5 ml of solution A to 99.5 ml of L.J culture media.

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
Preparation of Amikacin

Preparation of Amikacin solution (10mcg/ml)

∗ Dissolve 40 mg of amikacin in 10 ml of distilled water.
(AM concentration is 4000 mcg /ml= stock solution)
∗ Add 1 ml of stock –solution into 9 ml of distilled–water.
(AM concentration is 400 mcg/ml= solution A)
∗ Sterilize stock-solution by using 0.22 µl Millipore –filter.
∗ Add 1ml of filtrated solution, to 99 ml of L.J culture Media.
(AM concentration is 4 mcg/ml= working solution)

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Drug Containing Media

The drug containing medium is distributed in volumes of 6-8 ml in sterile 17mm × 170mm screw-capped test tubes, coagulated at 85^0C for 45 minutes ,and allowed to cool at room temperature for 24 hours; the screw caps are then tightened and the tubes stored at 4^0C. The period of validity of the media stored at 4^0C is 1 months.

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McFarland Standards
Barium Sulfate turbidity

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The McFarland Scale is a scale numbered from 1 to 10 which represents specific concentrations of bacteria/ml.

McFarland Scale	No. Bacteria ($\times 10^6$/ml)
1	300
2	600
3	900
4	1200
5	1500
6	1800
7	2100
8	2400
9	2700
10	3000

McFarland Standards are tubes labeled 1 through 10 and filled with suspensions of Barium salts. Each tube approximates the turbidity of bacterial solutions corresponding to the McFarland Scale number. Thus, tube 7 represents the turbidity of bacteria at a concentration of 2.1×10^{9}/ml.

If you have a culture and wish to determine its approximate population, its turbidity can be visually compared to a set of McFarland Standards. If its turbidity falls between tubes 7 and 8, then the number of bacteria/ml will be between 2.1 and 2.4 billion per ml. The advantage of these standards is that no incubation time or equipment needed to estimate bacterial numbers.

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Preparation of McFarland standard

0.5 McFarland Standard	**Sulfuric acid(H_2So_4) ,1% Barium Chloride($CL_2Ba,2H_20$) ,1%**	**995 ml 5ml**
1.0 McFarland Standard	**Sulfuric acid ,1% Barium Chloride,1%**	**990 ml 10ml**
2.0 McFarland Standard	**Sulfuric acid,1% Barium Chloride,1%**	**980ml 20ml**
3.0 McFarland Standard	**Sulfuric acid ,1% Barium Chloride,1%**	**970 ml 30ml**
4.0 McFarland Standard	**Sulfuric acid ,1% Barium Chloride,1%**	**960 ml 40 ml**

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
Put 8-10 glass beads in screw cap tube, add 200 µl of physiological serum or distilled waters, then autoclave the tube.

* Transfer 8-10 mycobacterium colonies to a sterile screw cap tube containing glass beads; homogenize it on mixer for 4-5 minutes. Take one ml of D.W and add the suspension into it gradually until the density adjust to that of a Mac Farland No.1

Prepare the suspension as follows:

1. Transfer 0.5 ml of suspension equal to Mac Farlands No.1, into tube contains 4.5 ml distilled water. In this way, 10^{-1} dilution has been prepared (solution B)
2. Transfer 0.1 ml of solution B to 9.9 distilled water (Solution C =10^{-3})
3. Transfer 0.1 ml of solution C to 9.9 distilled water (Solution D 10^{-5})

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Firstly, inoculate 200 µL of Mycobacterium suspension (with 10^{-1}, 10^{-3}, 10^{-5} concentration) on to each slant of control media and drug – containing media.

* Spread over the suspension on the whole surface of medium and incubated at 37^{0}C.

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Report the growth on control and drug containing media after 28, 42 days.

* First, read the mean of growth on control media and record it
* Record mean of growth of 10^{-1} tubes.
* Record mean of growth of 10^{-3} tubes.
* Record mean of growth of 10^{-5} tubes.

When no growth was observed on drug – containing media, the mycobacterium is sensitive.

* When enough growth was observed, count the number of colony, then calculates the percentage of resistance by this formula:

$$\%100 \times \frac{\text{Number of colonies on the drug}}{\text{Number of colonies on the control}}$$

Drug	Number of colonies		% resistant
	10^{-1}	10^{-3}	
Control	4+	60C	
INH 0.2 µg	120C	0	2%
SM 2. 0 µg	0	0	0
RIF 1.0µg	0	0	0

The 10^{-1} control inoculum gave 4+, or colonies too numerous to count. The hundredfold dilution (10^{-3}) gave 60 countable colonies. Therefore, the 4+ = 60 × 100 or 6000 colonies. Knowing this figure, the proportion of resistant bacilli may be calculated as 2% to 0.2 µg /ml INH.
If more than 1% of the test population is resistant to the drug under test, this suggests that resistance to that drug has developed.
An alternative way to calculate the percentage of resistant bacilli is:

$$\frac{\text{Colony count on drug medium at lower dilution } (10^{-1})}{\text{Colony count on control medium at higher dilution } (10^{-3})} = \text{\% resistance} \quad \frac{120}{60} = 2\%$$

Address: Mycobacteriology Research Centre, Iranian National Reference TB Laboratory, National Research Institute of Tuberculosis & Lung Disease (NRITLD), Shaheed Beheshti University of Medical Sciences, Shaheed Bahonar Ave, Darabad, 19556, P.O 19575/154, Tehran, Iran.
Tel: (0098-21) 20105555 or 20109505
Fax: (0098-21) 20109484 or 20109505
E-mail: mrc @ nritld.ac.ir, pfarnia @ hotmail.com, pfarnia @ nritld.ac.ir, Website: www.mrc.ac.ir

This document was created with Win2PDF available at http://www.daneprairie.com.
The unregistered version of Win2PDF is for evaluation or non-commercial use only.