In The Name Of God

Shaheed Beheshti University of
Medical sciences and Health Services

Collaborating Center for Tuberculosis
East Mediterranean Region

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National Research Institute
Tuberculosis and Lung Disease

Collaborating Center for Tuberculosis
Laboratory IUATLD

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Introduction

The Iranian Reference TB Laboratory receives 80 to 120 pulmonary and Extra-Pulmonary specimens per day.
The specimens were digested and decontaminated by petroff's (using 4% NAOH) method.
The procedures written in this manual based on procedures that performed in this laboratory.

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
Label a new and clean slide with the relevant code number

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Using the rough end of a broken bamboo stick, select and pick up the yellowish purulent particles of sputum.

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Spread the sputum on the glass slide evenly by smearing, repeatedly in coil-like patterns with approximately 1-2 cm in width and 2-3 cm in length.

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
Allow the smear to dry completely at room temperature (for at least 30 minutes)

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Allow the slide to fix on an incubator (65-75 ºC for 45 to 60 minute) or pass slide through burner flame as for other bacteriological smears.

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Label a new and clean slide with the relevant code number.

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After digestion and decontamination, prepare the smear from sediment. Spread the sputum on the glass slide evenly by smearing, repeatedly in coil-like patterns with approximately 1-2 cm in width and 2-3 cm in length.

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Allow the smear to dry completely at room temperature. (For at least 30 minutes)

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Allow the slide to fix on an incubator (65-75 ºC for 45 to 60 minute) or pass slide through burner flame as for other bacteriological examination.

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Discard the pasture pipette in disinfectant solution. (5%Phenol)

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✍ **Reagents**

Stock -Solution 1		Stock -Solution 2	
Fuchsin		**Phenol**	
Basic fuchsin	3.0g	Phenol Melted *	5 ml
Ethanol, 95%	100 ml	Distilled water	95ml

* (First melt phenol in 41 ºC incubator and then add it to distilled water)

Carbol fuchsin mix solution: 10 ml from Stock solution 1+ 90 ml of solution 2

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✍ **Decolourising agent: 3% acid-alcohol**

(Concentrated hydrochloric acid 3 ml, 95% ethanol 97 ml)

Carefully add 3ml of concentrated hydrochloric acid to 97 ml of 95% ethanol. Always add acid slowly to alcohol, not vice versa. Store the solution in bottle. Label bottle with name of reagent and dates of preparation the solution can be stored at room temperature for six to twelve months.

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✍ Counter stain: Methylene blue

Methylene blue chloride	0.3g
Distilled water	100 ml

Dissolve 0.3 g of methylene blue chloride in 100 ml of distilled water and store in an amber bottle. Label bottle with name of reagent and dates of preparation and expiry.

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
Place the slides on a staining rack (maximum 12). Ensure that slides do not touch each other.

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Flood entire slide with carbol fuchsin, which has been filtered prior to use.

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Heat the slide too slowly until it is steaming (Low or intermittent heat for three to five minutes); avoid boiling or drying of stain during heating.

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
Rinse each slide individually in a gentle stream of running water until all free stain is washed a way.

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Flood the slide with the decolourising solution (3% acid-alcohol) for maximum of three minutes.

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Rinse the slide thoroughly with water. Drain excess water from the slide.

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
Flood the slide with counter- stain for 60 seconds.

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

Rinse the slide thoroughly with water. Drain excess water from the slide; Allow smear to air dry.

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

Add one or two drops of immersion oil to the left edge of the stained smear.

If a smear is 2 × 3 cm in size
One (1) horizontal line is equivalent to 150 fields.
Two (2) vertical line equivalent to 100 fields.

If a smear is 1 × 2 cm in size
One (1) horizontal line is equivalent to 100 fields.
Two (2) vertical line equivalent to 100 fields.

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Reporting the results*

Recording

* No acid-fast bacilli found in at least **100** fields — **Negative**
* **1** to **9** acid-fast bacilli per **100** fields — **Exact number of bacilli**
* **10** to **99** acid-fast bacilli per **100** fields — +
* **1** to **10** acid-fast bacilli per field — ++
* **More** than **10** acid-fast bacilli per field — +++
* Consider basic morphology and color to accept a bacterium as a true acid-fast bacilli.
* If the smear preparations of fluid specimens are positive, prepare another smear from the same specimen and after careful examination report the final result.

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At the end of examination, take the slide from the microscopic stage then clean the objective lens with lens paper moistened with ethyl ether-alcohol. Petroleum benzene or xylene could be used instead of ethyl ether –alcohol.

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Keep the slides for quality control.

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If the result of smear microscopy is positive, record the code number, reader name, date and result on the label; stick the label on the slide; then store the slides in slide storage box in the sequence of code number.

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Auramine O mix solution: Mix solution 1 and 2, store in a tightly stopper amber bottle away from heat and light.

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Decolourising agent: 0.5 % acid-alcohol

(Concentrated hydrochloric acid 0.5 ml, 70% ethanol 100 ml)

Carefully add 0.5 ml of concentrated hydrochloric acid to 100 ml of 70% ethanol. Always add acid slowly to alcohol, not vice versa. Store the solution in bottle. Label bottle with name of reagent and dates of preparation the solution can be stored at room temperature for six to twelve months.

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✍ Counter stain: Potassium permanganate

Potassium permanganate	0.5g
Distilled water	100 ml

Dissolve 0.5 g of Potassium permanganate in 100 ml of distilled water and store in an amber bottle. Label bottle with name of reagent and dates of preparation and expiry. Store at room temperature for 3 months

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Flood entire slide with auramine O and allow to stain for 15 minutes,ensuring that stainig solution remains on smears. Do not heat and do not use filter paper strips

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Rinse the slide thoroughly with distilled water. Tap water contains chlorine which may interfere with fluorescence

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Decolourise with acid-alcohol for two minutes, Rinse with distilled water and drain.

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Flood smears with either potassium permanganate and allow counterstain for two minutes. Time is critical with potassium permanganate because counterstaining for longer may quench the fluorescence of acid –fast bacilli

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Allow smears to air – dry. Do not blot. Read as soon as possible after staining.

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Example: Suppose 20 acid - fast bacilli are observed per field using the 450X magnification. If this number is divided by the magnification factor of 4 according to table, the comparable number of bacilli that would have been observed under 1000X is 5 per field. The laboratory result should therefore read 2+

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The minimum number of acid-fast bacilli necessary to produce a positive smear result has been estimated to be between 5000 and 10000 per milliliter. The sensitivity can further be improved by examination of more than one smear from a patient. Because tuberculosis lesions in the lungs may drain intermittently, it is possible for a specimen to be negative on one day and positive the next; we collected **three specimens** for diagnosis.

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Removal of immersion oil from the slide
Prepare two jars of xylene
Dip the examined slide in the first jar for 30 minutes to remove the immersion oil and dip the slide again in the second jar to remove any remaining oil.
If xylene is not available, careful sweep using soft tissue paper can be done to remove oil from smear slide. Keep them in the slide box for confirmation or for quality check

Address: Mycobacteriology Research Centre, Iranian National Reference TB Laboratory, National Research Institute of Tuberculosis & Lung Disease (NRITLD), Shaheed Beheshti University of Medical Sciences, Shaheed Bahonar Ave, Darabad, 19556, P.O 19575/154, Tehran, Iran.
Tel: (0098-21) 20105555 or 20109505
Fax: (0098-21) 20109484 or 20109505
E-mail: mrc @ nritld.ac.ir, pfarnia @ hotmail.com, pfarnia @ nritld.ac.ir, Website: www.mrc.ac.ir

This document was created with Win2PDF available at http://www.daneprairie.com.
The unregistered version of Win2PDF is for evaluation or non-commercial use only.